

*The rapid swirling motion during vortex disrupts the liquid's surface, which is the primary point of aerosol creation. Also, vortexing generates turbulence that increases shear force at the liquid's surface.*

*The combination of turbulent flow and shear force separates particles from the liquid phase, launching them into the air as an aerosol.*

## VORTEX: Safety Practices

- Use vortex inside a biosafety cabinet or a fume hood whenever working with biohazardous materials. This confines any aerosols produced and provides a physical barrier against splash.
- Do not fill vessels to the brim, as this increases the chance of splashing and aerosol generation.
- Perform vortexing with sealed tubes or bottles whenever possible.
- Use vortex at the lowest speed necessary to mix the sample.



*Use of fluorescein to demonstrate material dispersion from an open tube mixed with a vortex.*

