
## Biosafety Office Broad Spectrum Surface Disinfectant – 10% Bleach



## 1 part (e.g. 100mL)

Ensure that your bleach has a concentration of **5 to 6%** of sodium hypochlorite.

> Bleach solutions are not very stable, and the antimicrobial activity of the chlorine typically decays over time. They should be **prepared fresh daily.**

Bleach should be stored in locations with temperatures between 50 and 70°F.



9 parts (e.g. 900mL)

To avoid degradation, store bleach solutions in containers that are opaque (to minimize light exposure) and closed (to minimize air exposure).